


Radio-Telemetry Monitoring of Blood Pressure in Conscious, Freely Moving, Undisturbed Mice


Key Features & Capabilities

- Applicable to all mouse strains and models
- Endpoints: Blood pressure and activity
- Continuous data collection for up to 3 months using the PA-C10 transmitter
- Capacity to monitor up to 60 individually housed mice
- Powerful DSI Dataquest™ A.R.T. software enables precise data acquisition and analysis
- Blood and tissue sampling is available

Diurnal Variation in Normal Mice



Anti-Hypertensive Effect in Mice



The Most Accurate and Precise System for Measuring Blood Pressure